


HUMAN AB SERUM: TRUE OFF-THE-CLOT **QUALITY TESTED SERUM FOR YOUR NEEDS**

**Human Male
AB Serum**
(Off-the-clot)

Akron's AB serum pools are produced from true off-the-clot (OTC) human serum, from healthy males who test negative for HBsAg, HCV and HIV 1/2 and is ideal for culturing cells associated with the human immune system.



AKRONBIOTECH

Human Male AB Serum (Off-the-clot)



FEATURES

- Male AB serum is low in estrogen and HCG
- Viral tested and found negative for HBsAg, HCV and HIV 1/2.
- Sterility testing is performed according to USP guidelines.
- Gamma irradiation and heat inactivation available upon request

Product	Cat#
Human AB Serum (True, off-the-clot)	AK9905
Human AB-Serum (Converted from plasma)	AK9340

RELATED PRODUCTS

- Platelet Lysate
- Human Serum
- Cord Blood Serum
- Affinity-stripped Plasma



Contact us at
sales@akronbiotech.com
for a sample or to receive
more information.



www.akronbiotech.com