


细胞培养产品

CELL CULTURE PRODUCTS

SCIENCE ORIENTED



上海道鹏生物科技有限公司
Shanghai XP Biomed Ltd.



C 目录 CONTENTS

VivaCell系列产品

牛血清	P02
其它动物血清	P04
完全培养基	P05
液体基础培养基	P06
液体缓冲液	P09
抗生素与消毒试剂	P11
细胞辅助溶液	P13
消化与冻存试剂	P14
无血清培养基与添加剂	P15
细胞与基因治疗培养基与添加物	P16
分子生物学试剂	P16
细胞遗传学	P17

Dr.Cell 希尔博士-技术服务

细胞健检 STR鉴定	P19
物种鉴定	
染色体核型分析	
其他细胞检测	
外泌体 外泌体NTA粒径、电位和浓度检测	
外泌体TEM检测	
外泌体WB检测	
外泌体提取	

PART ONE

VivaCell 系列产品 VivaCell Products



20年+
专注细胞生物学领域

VivaCell品牌由道鹏生物于2012年创立，涵盖：全线细胞培养血清和基础培养基系列、羊水细胞培养基、外周血淋巴细胞培养基、间充质干细胞分化染色试剂、组织消化液、细胞培养添加物、上皮 / 内皮细胞完全培养基、超低内毒素基础培养基 (SLE 系列) 等细胞培养产品。

01/ 牛血清



■ 甄选全球优质血源

血清均可追溯源头 | 来自无任何疯牛病疫情国家

■ 严苛工艺 多重质控

3次0.1µm过滤/6种以上细胞株培养测试 | 多种病毒/支原体/无菌检测

明星 产品



特级胎牛血清

适用于培养各种癌细胞株、胚胎细胞、白血病细胞、成纤维细胞、黑色素瘤细胞、巨噬细胞以及部分原代细胞。



精选胎牛血清

细胞培养增殖性能优越，适用多种细胞株培养；批间差小；多种病毒/支原体/无菌检测；低内毒素



胚胎干细胞专用胎牛血清

适用于培养诱导性多能干细胞、胚胎干细胞，人及小鼠干细胞、成体干细胞及原代细胞长期培养。



间充质干细胞专用胎牛血清

从多个批次的胎牛血清筛选出特定批次，经过骨髓，脂肪组织及脐带来源的间充质干细胞的细胞活性严格测试，通过传代及分化比例测试，可用于培养不同来源（脐带、脂肪、骨髓、牙髓或者胎盘等）的间充质干细胞，并可支持间充质干细胞后续分化成为硬骨、软骨和脂肪细胞等。

牛血清产品

产品名称 Product		规格 Size	货号 P/N	
胎牛血清	Certified Fetal Bovine Serum 特级胎牛血清	500ml 50ml 50ml*10支/盒	C04001-500 C04001-050 C04001-050X10	
	Fetal Bovine Serum 胎牛血清	500ml	C2910-0500	
	Certified Fetal Bovine Serum Qualified for Mesenchymal Cells 间充质干细胞专用胎牛血清	500ml 50ml 50ml*10支/盒	C04400-500 C04400-050 C04400-050X10	
	Certified Fetal Bovine Serum Qualified for Human Embryonic Stem Cell 胚胎干细胞专用胎牛血清	500ml 50ml 50ml*10支/盒	C04002-500 C04002-050 C04002-050X10	
	Fetal Bovine Serum-Select 精选胎牛血清	500ml 50ml 50ml*10支/盒	C2880-0500 C2880-0050 C2880-0050X10	
	Certified Gold Fetal Bovine Serum-Australia Origin 金牌胎牛血清-澳洲来源	500ml 50ml 50ml*10支/盒	C2820-0500 C2820-0050 C2820-0050X10	
	Certified Gold Fetal Bovine Serum-New Zealand Origin 金牌胎牛血清-新西兰来源	500ml 50ml 50ml*10支/盒	C2830-0500 C2830-0050 C2830-0050X10	
	其他牛血清	Adult Bovine Serum 成年牛血清	500ml	C2610-0500
		Newborn Calf Serum 新生牛血清	500ml 100ml	C2620-0500 C2620-0100
		Neonatal Calf Serum 类胎牛血清	500ml 100ml	C2630-0500 C2630-0100
		特殊加工工艺血清	Certified Fetal Bovine Serum Qualified for Tetracycline Regulated Systems 四环素调控系统胎牛血清	500ml 50ml 50ml*10支/盒
	Certified Fetal Bovine Serum Qualified for Low IgG 低IgG胎牛血清		500ml 50ml 50ml*10支/盒	C2730-0500 C2730-0050 C2730-0050X10
Exosome Depleted Fetal Bovine Serum 去外泌体胎牛血清	500ml 50ml 50ml*10支/盒		C3801-0100 C3801-0050 C3801-0050X10	
Dialyzed Fetal Bovine Serum 透析胎牛血清	500ml 50ml 50ml*10支/盒		C3820-0500 C3820-0050 C3820-0050X10	
Charcoal Stripped Fetal Bovine Serum 碳吸附胎牛血清	500ml 50ml 50ml*10支/盒		C3830-0500 C3830-0050 C3830-0050X10	

产品名称 Product

规格 Size

货号 P/N

特殊加工 工艺血清	Lipid Depleted Fetal Bovine Serum 去脂胎牛血清	500ml 50ml 50ml*10支/盒	C3840-0500 C3840-0050 C3840-0050X10
	Certified Fetal Bovine Serum Heat Inactivated 热灭活特级胎牛血清	500ml 50ml	C04001-500HI C04001-50HI
	Certified Fetal Bovine Serum Qualified for Human Embryonic Stem Cell Heat Inactivated 热灭活胚胎干细胞专用胎牛血清	500ml	C04002-500HI
	Certified Fetal Bovine Serum Qualified for Tetracycline Regulated Systems Heat Inactivated 热灭活四环素调控系统胎牛血清	500ml	C2720-0500HI
	Newborn Calf Serum Heat Inactivated 热灭活新生牛血清	500ml	C2620-0500HI

其它动物血清产品

产品名称 Product

规格 Size

货号 P/N

其他 动物血清	Porcine Serum 猪血清	1000ml 500ml	C2500-1000 C2500-0500
	Horse Serum 马血清	500ml 100ml	C2510-0500 C2510-0100
	Rabbit Serum 兔血清	100ml	C2520-0100
	Goat Serum 羊血清	100ml	C2530-0100
	Donkey Serum 驴血清	500ml 100ml	C2540-0500 C2540-0100
	Chicken Serum 鸡血清	100ml	C2550-0100
	Mouse Serum 小鼠血清	100ml 50ml	C2560-0100 C2560-0050
	Rat Serum 大鼠血清	100ml 50ml	C2570-0100 C2570-0050
	Horse Serum Heat Inactivated 热灭活马血清	500ml	C2510-0500HI

备注:全线血清产品可以提供伽马照射服务与热灭活处理

★ More products to come !!!

02/ 完全培养基 液体基础培养基 液体缓冲液



■ 完全培养基

- 营养全面，无需额外添加成分
- 即用型，实现高效增殖，为细胞生长提供稳定微环境

■ 液体基础培养基

- 种类丰富 | 支持个性化定制 | 一站式满足多种细胞培养
- 标准化基础营养配方，为细胞生长提供核心支撑
- 多品类覆盖贴壁/悬浮、原代/细胞系，适配不同细胞代谢需求
- 覆盖科研实验、规模化生产等多场景，适用细胞培养全流程

■ 液体缓冲液

- 无需溶解，即用型，减少操作误差
- 精准pH缓冲范围，为实验体系提供稳定酸碱环境

完全培养基产品

产品名称 Product		规格 Size	货号 P/N
完全培养基	Endothelial Cell Complete Medium	500ml	C3660-0500
	内皮细胞完全培养基	100ml	C3660-0100
	Epithelial Cell Complete Medium	500ml	C3650-0500
	上皮细胞完全培养基	100ml	C3650-0100
	Expander Complete Single Cell Culture Medium	500ml	C3670-0500
	单细胞增殖完全培养基		

液体基础培养基产品

产品名称 Product	规格 Size	货号 P/N
RPMI 1640	RPMI 1640, with L-Alanyl-L-Glutamine RPMI 1640 (含 L-丙氨酰-L-谷氨酰胺)	500ml C3001-0500
	RPMI 1640, with L-Glutamine RPMI1640 (含谷氨酰胺)	500ml C3010-0500
	RPMI 1640, with L-Glutamine, SLE Grade 超低内毒素RPMI1640 (含谷氨酰胺), SLE级别	500ml C3011-0500
	RPMI Medium 1640, w/o D-Glucose, w/o L-Glutamine RPMI1640 (不含葡萄糖, 不含谷氨酰胺)	500ml C3015-0500
	RPMI Medium 1640, w/o Phenol Red, w/o L-Glutamine RPMI1640 (不含酚红, 不含谷氨酰胺)	500ml C3016-0500
	RPMI Medium 1640, w/o L-Glutamine RPMI1640 (不含谷氨酰胺)	500ml C3017-0500
	RPMI Medium 1640, with 25mM Hepes, with L-Glutamine RPMI1640 (含Hepes, 含谷氨酰胺)	500ml C3018-0500
	RPMI Medium 1640, w/o Phenol Red, with L-Glutamine RPMI1640 (不含酚红, 含谷氨酰胺)	500ml C3019-0500
MEM-E	Minimum Essential Medium Eagle,(MEM-E), with Earle's salts and L-Glutamine, w/o NEAA MEM-E (含谷氨酰胺)	500ml C3050-0500
	Minimum Essential Medium (MEM), Earle's Salts Base, with Non-Essential Amino Acids, w/o L-Glutamine MEM (含非必需氨基酸NEAA, 不含谷氨酰胺)	500ml C3053-0500
	Minimum Essential Medium (MEM), Earle's Salts Base, with Non-Essential Amino Acids, w/o L-Glutamine MEM (含非必需氨基酸NEAA, 含谷氨酰胺)	500ml C3054-0500
MEM-α	Minimum Essential Medium Alpha (MEM-α), with L-Glutamine, w/o Ribonucleosides, Deoxyribonucleosides MEM-α (含谷氨酰胺, 不含核苷核酸和脱氧核苷核酸)	500ml C3060-0500
	Minimum Essential Medium Alpha (MEM-α), with L-Glutamine, w/o Ribonucleosides, Deoxyribonucleosides, SLE Grade 超低内毒素MEM-α (含谷氨酰胺, 不含核苷核酸和脱氧核苷核酸), SLE级别	500ml C3061-0500
	Minimum Essential Medium Alpha (MEM-α), with L-Glutamine, w/o Ribonucleosides, Deoxyribonucleosides, Phenol red MEM-α (含谷氨酰胺, 不含核苷核酸和脱氧核苷核酸, 酚红)	500ml C3062-0500
	Minimum Essential Medium Alpha (MEM-α), with L-Glutamine, w/o Ribonucleosides, Deoxyribonucleosides, w/o Phenol red, SLE Grade 超低内毒素MEM-α (含谷氨酰胺, 不含酚红), SLE级别	500ml C3063-0500
	Minimum Essential Medium Alpha (MEM-α), w/o Ascorbic acid, with L-Glutamine, with Ribonucleosides, Deoxyribonucleosides MEM-α (不含抗坏血酸, 含谷氨酰胺, 含核苷核酸和脱氧核苷核酸)	500ml C3064-0500

产品名称 Product
规格 Size 货号 P/N

MEM-α	Minimum Essential Medium Alpha (MEM-α), w/o Ascorbic acid, Ribonucleosides, Deoxyribonucleosides, with L-Alanyl-Glutamine, Phenol red	500ml	C3065-0500
	MEM-α (不含抗坏血酸, 核苷核酸和脱氧核苷核酸, 含L-丙氨酰-L-谷氨酰胺, 酚红)		
	Minimum Essential Medium Alpha (MEM-α), w/o Ascorbic acid, w/Ribonucleosides, Deoxyribonucleosides, L-Alanyl-Glutamine, Phenol red, SLE Grade	500ml	C3002-0500
	超低内毒素 MEM-α (不含抗坏血酸, 含核苷核酸和脱氧核苷核酸, L-丙氨酰-L-谷氨酰胺, 酚红) , SLE级别		
	Minimum Essential Medium-Alpha (MEM-α), w/o Ascorbic acid, Ribonucleosides, Deoxyribonucleosides ,with L-Glutamine, Phenol red, SLE Grade	500ml	C3066-0500
	超低内毒素MEM-α (不含抗坏血酸, 核苷核酸和脱氧核苷核酸, 含谷氨酰胺, 酚红) , SLE级别		
	Minimum Essential Medium-Alpha (MEM-α), w/o Ascorbic acid, Ribonucleosides , Deoxyribonucleosides , Phenol red , with L-Glutamine, SLE Grade	500ml	C3067-0500
	MEM-α (不含抗坏血酸, 核苷核酸和脱氧核苷核酸, 酚红, 含谷氨酰胺)		
	Minimum Essential Medium-Alpha (MEM-α), w/o Ascorbic acid, Ribonucleosides, Deoxyribonucleosides, with L-Glutamine, Phenol red	500ml	C3068-0500
	MEM-α (不含抗坏血酸, 核苷核酸和脱氧核苷核酸, 含谷氨酰胺, 酚红)		
	Minimum Essential Medium-Alpha (MEM-α), w/o Ascorbic acid, Ribonucleosides, Deoxyribonucleosides, Phenol red, with L-Glutamine	500ml	C3069-0500
	MEM-α (不含抗坏血酸, 核苷核酸和脱氧核苷核酸, 酚红, 含谷氨酰胺)		
L-15	Leibovitz's L-15 Medium, with L-Glutamine	500ml	C3070-0500
	L-15 (含谷氨酰胺)		
	Leibovitz's L-15 Medium with L-Glutamine, w/o Sodium bicarbonate, Phenol Red	500ml	C3071-0500
	L-15 (含谷氨酰胺, 不含碳酸氢钠, 酚红)		
DMEM高糖	Dulbecco' s Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Alanyl-L-Glutamine, w/o Sodium Pyruvate	500ml	C3101-0500
	DMEM高糖 (含 L-丙氨酰-L-谷氨酰胺, 不含丙酮酸钠)		
	Dulbecco's Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Alanyl-L-Glutamine, with Sodium Pyruvate	500ml	C3103-0500
	DMEM高糖 (含 L-丙氨酰-L-谷氨酰胺, 含丙酮酸钠)		
	Dulbecco' s Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Glutamine, w/o Sodium Pyruvate	500ml	C3110-0500
	DMEM高糖 (含谷氨酰胺, 不含丙酮酸钠)		
	Dulbecco' s Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Glutamine, w/o Sodium Pyruvate, SLE Grade	500ml	C3111-0500
	超低内毒素DMEM高糖 (含谷氨酰胺, 不含丙酮酸钠), SLE级别		

产品名称 Product
规格 Size 货号 P/N

DMEM高糖	Dulbecco's Modified Eagle Medium (DMEM), With 4.5g/L D-Glucose (High Glucose), with L-Glutamine, w/o Sodium Pyruvate, w/o Phenol Red DMEM高糖 (含谷氨酰胺, 不含丙酮酸钠, 不含酚红)	500ml	C3112-0500
	Dulbecco's Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Glutamine, with Sodium Pyruvate DMEM高糖 (含谷氨酰胺, 含丙酮酸钠)	500ml	C3113-0500
	Dulbecco's Modified Eagle Medium (DMEM), with 4.5g/L D-Glucose (High Glucose), with L-Glutamine, with Sodium Pyruvate, SLE Grade 超低内毒素DMEM高糖 (含谷氨酰胺, 含丙酮酸钠, SLE级别)	500ml	C3114-0500
	Dulbecco's Modified Eagle Medium (DMEM), High Glucose, w/o Sodium Pyruvate, w/o L-Glutamine DMEM高糖 (不含丙酮酸钠, 不含谷氨酰胺)	500ml	C3115-0500
	Dulbecco Modified Eagle Medium (DMEM) 2X Conc., High Glucose, w/o L-Glutamine, w/o Sodium Pyruvate, w/Sodium Bicarbonate DMEM高糖 (不含谷氨酰胺和丙酮酸钠, 含碳酸氢钠, 2倍浓缩液)	500ml	C3117-0500
	Dulbecco Modified Eagle Medium (DMEM), High Glucose, w/o L-Glutamine, Sodium Pyruvate, Phenol red DMEM高糖 (不含谷氨酰胺和丙酮酸钠, 酚红)	500ml	C3118-0500
DMEM低糖	Dulbecco's Modified Eagle Medium (DMEM), with 1g/L D-Glucose (Low Glucose), with L-Glutamine, with Sodium Pyruvate DMEM低糖 (含谷氨酰胺)	500ml	C3120-0500
DMEM无糖	Dulbecco's Modified Eagle Medium (DMEM), w/o D-Glucose, w/o Sodium Pyruvate, w/o L-Glutamine DMEM (不含葡萄糖, 不含丙酮酸钠, 不含谷氨酰胺)	500ml	C3122-0500
DMEM: F12	Dulbecco's Modified Eagle Medium (DMEM): Nutrient Mixture F-12(Ham's) (1:1) , with L-Glutamine and Sodium Pyruvate, with 15mM HEPES DMEM: F12 (1:1) (含谷氨酰胺, 含丙酮酸钠, 含HEPES)	500ml	C3130-0500
	Dulbecco's Modified Eagle Medium (DMEM): Nutrient Mixture F-12 (Ham's) (1:1) with L-Glutamine and Sodium Pyruvate, with 15mM HEPES, SLE Grade 超低内毒素DMEM: F12 (1:1) (含谷氨酰胺, 含丙酮酸钠, 含HEPES), SLE级别	500ml	C3131-0500
	Dulbecco's Modified Eagle Medium (DMEM): Nutrient Mixture F-12 (Ham's) (1:1) with L-Glutamine, sodium Pyruvate, w/o HEPES DMEM: F12 (含谷氨酰胺, 含丙酮酸钠, 不含HEPES)	500ml	C3132-0500
	Dulbecco's Modified Eagle Medium (DMEM): Nutrient Mixture F-12 (Ham's) (1:1) with L-Glutamine, sodium Pyruvate, 15mM HEPES, w/o Pheno red DMEM: F12 (含谷氨酰胺, 含丙酮酸钠, 含 HEPES, 不含酚红)	500ml	C3133-0500
昆虫基础培养基	Grace's Insect Medium 昆虫培养基	500ml	C3150-0500
	Schneider's Drosophila Medium, with L-Glutamine 施耐德的果蝇培养基 (含谷氨酰胺)	500ml	C3160-0500

产品名称 Product

规格 Size

货号 P/N

其他基础培养基	Nutrient Mixture F-12 (Ham's), with L-Glutamine F-12 (Ham's) (含谷氨酰胺)	500ml	C3134-0500
	Nutrient Mixture F-12 (Ham's), w/o L-Glutamine F-12 (Ham's) (不含谷氨酰胺)	500ml	C3136-0500
	MCDB-153 (Modified), with L-Glutamine MCDB-153 (含谷氨酰胺)	500ml	C3170-0500
	Waymouth's MB 752/1 Medium Waymouth's培养基	500ml	C3175-0500
	TransfMEM- α Reduced Serum Medium 减血清细胞转染培养基	500ml	C3310-0500
	Glasgow's MEM (GMEM)	500ml	C3056-0500
	McCoy's 5A Medium (Modified), with L-Glutamine McCoy's 5A (含谷氨酰胺)	500ml	C3020-0500
	Medium M-199 (M199), Earle's Salts Base, with L-Glutamine M-199 (含谷氨酰胺)	500ml	C3030-0500
	Iscove's Modified Dulbecco's Medium (IMDM), with L-Glutamine, w/o Alpha-Thioglycerol, Beta-Mercaptoethanol IMDM (含谷氨酰胺)	500ml	C3040-0500
	Iscove's Modified Dulbecco's Medium (IMDM), with L-Glutamine, w/o Alpha-Thioglycerol, w/o Beta-Mercaptoethanol, w/o Phenol Red IMDM (含谷氨酰胺, 不含酚红)	500ml	C3042-0500

液体缓冲液产品

产品名称 Product

规格 Size

货号 P/N

HBSS	Hanks' Balanced Salt Solution HBSS平衡盐溶液 (含钙, 镁和酚红)	500ml	C3570-0500
	Hanks' Balanced Salt Solution, w/o Calcium and Magnesium, w/o Phenol Red HBSS平衡盐溶液 (无钙, 镁和酚红)	500ml	C3571-0500
	Hanks' Balanced Salt Solution, with Calcium and Magnesium without Phenol Red HBSS平衡盐溶液 (含钙和镁, 无酚红)	500ml	C3572-0500
	Hanks' Balanced Salt Solution, without Calcium and Magnesium, with Phenol Red HBSS平衡盐溶液 (无钙和镁, 有酚红)	500ml	C3573-0500
	Phosphate Buffered Saline (PBS), w/o Calcium and Magnesium PBS (不含钙, 镁)	500ml	C3580-0500
PBS	Phosphate Buffered Saline (PBS), w/o Calcium and Magnesium, SLE Grade 超低内毒素PBS (不含钙, 镁), SLE级别	500ml	C3581-0500
	Phosphate Buffered Saline (PBS), 10X Conc., w/o Calcium and Magnesium PBS缓冲液 (不含钙, 镁, 10倍浓缩液)	500ml	C3582-0500

产品名称 Product
规格 Size
货号 P/N

DPBS	Dulbecco's Phosphate Buffered Saline (DPBS), w/o Calcium and Magnesium DPBS (不含钙, 镁)	500ml	C3590-0500
	Dulbecco's Phosphate Buffered Saline (DPBS), w/o Calcium and Magnesium, SLE Grade 超低内毒素DPBS (不含钙, 镁), SLE级别	500ml	C3591-0500
	Dulbecco's Phosphate Buffered Saline (DPBS), 10X Conc., w/o Calcium and Magnesium DPBS缓冲液 (不含钙, 镁和酚红, 10倍浓缩液)	500ml	C3592-0500
	Dulbecco's Phosphate Buffered Saline (DPBS), with Calcium and Magnesium, w/o Phenol Red DPBS缓冲液 (含钙和镁, 不含酚红)	500ml 100ml	C3593-0500 C3593-0100
其它	HEPES Buffer Solution, pH 7.3 at 37°C, 1.0 Molar HEPES缓冲液	100ml	C3544-0100
	Earle's Balanced Salt Solution EBSS缓冲液	500ml	C3548-0500
	Earle's Balanced Salt Solution, with Calcium and Magnesium with Phenol Red EBSS缓冲液 (含钙和镁, 含酚红)	500ml	C3549-0500
	Alsever's solution Alsever 溶液	500ml	C3556-0500
	Ethylenediaminetetraacetic Acid (EDTA) Disodium Salt Solution, 0.05% in DPBS EDTA溶液	100ml	C3550-0100
	Water, Cell Culture Grade, SLE Grade 水 (细胞培养级别), SLE级别	500ml	C3560-0500
	Highly Purified Water (DNase and RNase-free), SLE Grade, NTA Passed 超纯水, SLE级别, NTA检测通过	500ml	C3562-0500

★ More products to come !!!

03/ 抗生素与消毒试剂



细胞健康守护卫士 | 卓越品质 护航您的科研成果

细胞培养中使用的抗生素包括单抗，双抗，三抗以及支原体清除试剂，可用于预防或控制微生物污染，确保细胞的健康生长和实验的可靠性。

明星 产品



青链双抗（货号：C3420）

- ★ 广谱杀菌活性
- ★ 培养细胞最常用的抗生素



两性霉素B（货号：C3430）

- ★ 抗真菌谱范围广
- ★ 抗真菌活性高



培养箱水盘消毒剂（货号：C3480）

- ★ 广谱抗菌活性
- ★ 稀释后即，对不锈钢无腐蚀性

特别 推荐



灭菌喷液（货号：C3490）

- ★ 可预防真菌及其孢子，细菌及其孢子，支原体和病毒（包括HIV和HBV）
- ★ 可生物降解、无毒、无刺激，对设备无腐蚀



支原体清除试剂系列

- ★ 无菌过滤(0.1 μm)
- ★ 清除支原体污染无死角



抗生素与消毒试剂产品

抗生素系列

产品名称 Product	规格 Size	货号 P/N	抗菌谱		
			G+细菌	G-细菌	真菌
青链双抗 Penicillin-Streptomycin Solution	100ml	C3420-0100	√	√	
青链双抗, 10X浓缩液 Penicillin-Streptomycin Solution, 10X Conc.	100ml	C3421-0100	√	√	
青链霉素, 制霉菌素三抗溶液 Penicillin-Streptomycin Nystatin Solution	100ml	C3422-0100	√	√	√
青链霉素, 两性霉素B三抗 Penicillin-Streptomycin Amphotericin B Solution	100ml	C3423-0100	√	√	√
青链霉素, 新霉素三抗 Penicillin-Streptomycin Neomycin Solution	100ml	C3424-0100	√	√	
两性霉素B Amphotericin B Solution	100ml	C3430-0100			√
制霉菌素混悬液 Nystatin Suspension	20ml	C3440-0020			√
硫酸庆大霉素溶液 Gentamycin Sulfate Solution	100ml	C3450-0100			√
卡那霉素溶液 Kanamycin Sulphate Solution	100ml	C3460-0100			√

支原体清除试剂

产品名称 Product	规格 Size	货号 P/N
支原体处理试剂1 MycoNator™-1 Antibiotic Solution 100X Conc.	20ml	C3470-0020
支原体处理试剂2 MycoNator™-2 Antibiotic Solution 100X Conc.	20ml	C3471-0020
支原体处理试剂3 MycoNator™-3 Antibiotic Solution 100X Conc.	100ml	C3472-0100
	20ml	C3472-0020
	10ml	C3472-0010

设备仪器预防试剂/其它

产品名称 Product	规格 Size	货号 P/N
培养箱水盘消毒剂 Akouguard-1, Water tray in CO ₂ incubators Disinfection Solution	100ml	C3480-0100
水浴锅消毒剂 Akouguard-2, Water bath Disinfecting Solution	500ml	C3481-0500
	50ml	C3481-0050
灭菌喷液 Pharmacidal, for disinfecting surfaces	1L	C3490-1000
	100ml	C3490-0100
MycoEZ™ 支原体qPCR检测试剂盒 MycoEZ™, Mycoplasma qPCR Kit	25rxn	M5010-0025

04/细胞辅助溶液



间充质干细胞成骨/成脂/成软骨分化试剂盒是一款用于间充质干细胞诱导成骨分化的产品，可用于骨髓、脐带、脂肪等组织来源的间充质干细胞的诱导成骨/成脂/成软骨分化。台盼蓝溶液可用于检测细胞的活性，广泛应用于细胞培养、药物毒性测试、细胞移植等实验研究。

细胞辅助溶液产品

产品名称 Product	规格 Size	货号 P/N
Adipo-dye for adipocytes staining kit in adipogenic differentiation 间充质干细胞成脂分化染色试剂盒	150ml	C37A0-0150
Chondro-dye for chondrocytes staining kit in chondrogenic differentiation 间充质干细胞成软骨分化染色试剂盒	150ml	C37B0-0150
Osteo-dye for osteocytes staining kit in osteogenic differentiation 间充质干细胞成骨分化染色试剂盒	150ml	C37C0-0150
Trypan Blue Solution, 4mg/ml in Saline 台盼蓝溶液	100ml	C3730-0100

05/ 消化与冻存试剂



胰酶用于解离细胞；组织消化液能分解组织获取单细胞悬液；冻存液用于细胞冻存，助力细胞实验及培养流程。

明星产品



重组胰酶（含EDTA）（货号：C3536）

Recombinant Trypsin EDTA Solution

- 批次差异小，纯度高
- 无动物源，无病毒与支原体的污染风险



细胞冻存液3号（货号：C3525/C3526）

Cell Cryopreservation Solution 3

- 批次差异小
- 低内毒素，无污染
- 细胞存活率高

产品名称 Product

规格 Size 货号 P/N

组织消化液	Tissue Digestive Mix-I	100ml	C3502-0100
	组织消化液（RUO）	20ml	C3502-0020
	Cell Dissociation Solution (non-enzymatic)	100ml	C3505-0100
细胞冻存液	细胞分离溶液(无酶法)		
	Cell Cryopreservation Solution I	500ml	C3510-0500
	细胞冻存液 I号（无动物源成分）	100ml	C3510-0100
	Cell Cryopreservation Solution I, w/o phenol-red	500ml	C3511-0500
	细胞冻存液 I号（无动物源成分，无酚红）	100ml	C3511-0100
	Cell Cryopreservation Solution 3	100ml	C3525-0100
胰酶	无血清细胞冻存液 3号（无酚红）	50ml	C3525-0050
	Cell Cryopreservation Solution 3	100ml	C3526-0100
	无血清细胞冻存液 3号（有酚红）	50ml	C3526-0050
	Trypsin EDTA Solution A (0.25%), EDTA (0.02%), with Phenol Red	500ml	C3530-0500
	胰酶（胰蛋白酶0.25%，EDTA0.02%）	100ml	C3530-0100
	Trypsin EDTA Solution B (0.25%), EDTA (0.05%), with Phenol Red	500ml	C3532-0500
	胰酶（胰蛋白酶0.25%，EDTA0.05%）	100ml	C3532-0100
	Trypsin EDTA Solution C (0.05%), EDTA (0.02%), with Phenol Red	500ml	C3533-0500
	胰酶（胰蛋白酶0.05%，EDTA0.02%）	100ml	C3533-0100
	Trypsin Solution B (0.25% in DPBS), without Calcium, Magnesium, without Red Phenol	500ml	C3538-0500
	胰蛋白酶（不含钙和镁，不含酚红，不含EDTA）	100ml	C3538-0100
	Recombinant Trypsin EDTA Solution	500ml	C3536-0500
	重组胰酶（含EDTA）	100ml	C3536-0100
	Soybean Trypsin Inhibitor 50X Conc, 5mg/ml	20ml	C3553-0020
	大豆胰蛋白酶抑制剂		

06/ 无血清培养基与添加剂



额外添加细胞培养过程中所需的氨基酸、谷氨酰胺、丙酮酸钠、人胰岛素（重组）、人转铁蛋白（重组）、人白蛋白（重组）等。

无血清培养基与添加剂产品

产品名称 Product	规格 Size	货号 P/N
ITS Animal Component free (ACF) Cell Culture Medium Supplement (100X) ITS细胞添加物 (100X)	5g	C3210-0005
ITS-E Cell Animal Component Free (ACF) Culture Medium Supplement (100X) ITS-E细胞添加物 (100X)	5g	C3213-0005
MEM Non-Essential Amino Acids Solution, 100X Conc. MEM非必需氨基酸100X 浓缩	100ml	C3240-0100
MEM Amino Acids Solution, 50X Conc., w/o L-Glutamine MEM氨基酸50X 浓缩	100ml	C3250-0100
L-Glutamine Solution (200 mM) L-谷氨酰胺	500ml 100ml	C3543-0500 C3543-0100
Sodium Pyruvate Solution, 11.0mg/ml (100 mM) 丙酮酸钠	100ml	C3546-0100
Stable Glutamine (L-Alanyl-L-Glutamine) 长效型谷氨酰胺 (L-丙氨酰-L-谷氨酰胺)	500ml 100ml	C3547-0500 C3547-0100
EaSulin™ Recombinant Human Insulin, Powder 重组人胰岛素 (干粉)	1g 500mg 100mg	C6010-1000 C6010-0500 C6010-0100
EastFerrin™ Recombinant Human Transferrin, Powder 重组人转铁蛋白 (干粉)	5000mg 1000mg 500mg 100mg	C6510-5000 C6510-1000 C6510-0500 C6510-0100
EastBumin™ Recombinant Human Albumin, Powder 重组人白蛋白 (干粉)	100mg 5g 1g	C6520-0100 C6520-0005 C6520-0001

07/ 细胞与基因治疗 培养基与添加剂 分子生物学试剂



可有效支持干细胞等人类细胞生长。含有各类营养成分、生长因子与特殊添加物，确保细胞生长、增殖和功能表达。

细胞与基因治疗培养基与添加剂产品

产品名称 Product		规格 Size	货号 P/N
血小板裂解物	VivaPL™ Human Platelet Lysate , Heparin Free (RUO) 血小板裂解物 (科研级)	500ml	C8010-0500
		100ml	C8010-0100
		10ml	C8010-0010
	VivaPL™ Human Platelet Lysate, Heparin needed (GMP grade) 血小板裂解物, 需抗凝剂 (GMP级)	500ml	C8013-0500
		100ml	C8013-0100
		10ml	C8013-0010
脐带来源间充质干细胞外泌体	Exometics® UC-MSC Exosome 脐带来源间充质干细胞 (UC-MSC)外泌体, 浓度(particles/ml): 1*1010/ml	100ml	C6610-0100
		10ml	C6610-0010

一款用于快速且均匀将细菌涂布在培养基表面的玻璃珠，有利于细菌的均匀生长。

分子生物学试剂产品

产品名称 Product	规格 Size	货号 P/N
涂布玻璃珠 (Glass Plating Beads)	150g	M5800-0150

08/ 细胞遗传学



细胞遗传学产品用于培养骨髓、羊水、外周血来源的细胞，细胞增殖后可检测细胞染色体数目、结构异常及相关基因位点变化，能辅助疾病诊断、遗传咨询与科研。

明星产品

细胞培养基（货号：C3620/C3621）AmnioType-1 Medium / AmnioType-2 Medium

羊水培养基AmnioType是专门针对羊水和绒毛细胞培养而研发的，用于胎儿染色体核型分析。

- ①专门针对羊水和绒毛细胞培养设计
- ②完全培养基，缓冲效果更佳，适用于封闭和开放培养系统

骨髓细胞培养基（货号：C3630）Bone Marrow Karyotyping Medium

用于骨髓体外增殖培养，不具备对细胞的选择、诱导、分化功能。培养后用于临床诊断。骨髓培养基是一种完全补充的培养基，用于骨髓和其他造血细胞的短期培养，用于细胞遗传学研究和体外诊断。

- ①稳定的培养效果
- ②良好的细胞形态，显带效果好，染色体清晰，体积大，形态好，适用性强

细胞培养基（货号：C3601）PBType-1 Cell Culture Medium

用于外周血淋巴细胞体外增殖培养，不具备对细胞的选择、诱导、分化功能。使用该培养基培养增殖后，提供后续检测，如染色体核型分析。

- ①专利丁基胶塞设计-加样方便，收集和染色体制备在同一管中进行，无需开盖，避免污染
- ②HEPES缓冲体系-适用于封闭式细胞培养兼容自动化操作-适用于HANABI、Chromprep II等多种自动化核型分析仪

细胞同步化试剂盒（货号：C3611）Cell Syncornization Kit

细胞同步化培养技术是利用药物终止DNA合成，使细胞周期停止在S期。以便于稳定获取大量高质量中期分裂相。从而同步进入细胞分裂期，此时可以添加秋水仙碱，达到细胞周期同步化的目的。

- ①毒性低，处理时间短
- ②操作方便，培养时间灵活

细胞遗传学产品

产品名称 Product	规格 Size	货号 P/N
细胞培养基 AmnioType-1 Medium	100ml	C3620-0100
细胞培养基 AmnioType-2 Medium	100ml	C3621-0100
骨髓细胞培养基 MarrowType Bone Marrow Karyotyping Medium	100ml	C3630-0100
细胞培养基 PBType-1 Cell Culture Medium	5ml	C3601-0005
	100ml	C3601-0100
	5ml*30	C3601-0150
	500ml	C3601-0500
秋水酰胺溶液 Colcemid (Demecolcine) Solution, 10 ug/mL in DPBS	10ml	C3541-0010
植物血球凝集素 (液体) Phytohemagglutinin-M (PHA-M), Liquid	5ml	C3545-0005
氯化钾低渗处理液 Potassium Chloride Solution 0.075 M	100ml	C3540-0100
	500ml	C3540-0500
胰酶消化溶液 Trypsin EDTA Solution	100ml	C3680-0100
	500ml	C3680-0500
胰酶分带溶液 Trypsin Solution B	100ml	C3683-0100
	500ml	C3683-0500
吉姆萨染色液 Giemsa Stain	100ml	C3720-0100
	500ml	C3720-0500
细胞同步化试剂盒 Cell Syncornization Kit	8.8ml	C3611-0009

PART TWO



Dr.Cell 希尔博士-技术服务 Dr.Cell Technical Service

■ 外泌体NTA 粒径、电位和浓度检测

国内首家外泌体NTA检测服务，最快当天出报告

- 每个样本3重复检测，符合统计学意义
- 配备德国PMX纳米颗粒跟踪，分析仪Zetaview
- 样本处理简单，保持样本原始状态
- 可检测带有荧光抗体，标记的样本
- 提供外泌体粒径和绝对浓度数据
- 可检测zeta电位



■ TEM 电镜检测

作为外泌体鉴定的“金标准”之一，TEM 检测可明确样品中是否存在具有典型外泌体形态特征（如杯状或茶托状、具有双层膜结构）的囊泡，为外泌体的存在提供直接的形态学证据。

- 清晰呈现外泌体的完整、圆形或杯状形态，排除非囊泡结构（如蛋白质聚集体）的干扰
- 通过观察视野内是否存在明显的细胞碎片或杂质，辅助判断样品的纯度
- 提供具有发表级质量的高清TEM图像，可直接用于学术论文、报告及展示

■ 蛋白标志物检测（WB）

外泌体特异性标志蛋白检测服务，是为科研或临床转化需求提供的专业化技术服务，核心是通过标准化实验流程，精准检测外泌体标志性蛋白，验证其身份、纯度及表达差异。

- 结果可靠性与重复性更强，提供5种抗体指标（CD9/CD63/CD81/TSG101/Calnexin）
- 提供专业解读结果与报告

■STR细胞鉴定

短串联重复序列（STR）鉴定是基于人类基因组中存在的大量短串联重复序列进行的身份识别技术。不同个体的STR位点在重复单元的数量上存在显著差异，这种差异具有高度的个体特异性和遗传稳定性。

- 检测经验丰富：检测样本数量1000+，拥有丰富的细胞鉴定经验
- 辨识度高：人源细胞鉴定包含21个基因座位点，具有更高辨识度
- 准确性高：检测数据与国际权威数据库进行比对，以确保鉴定结果的准确性

■核型分析-人源细胞

iPSC核型分析服务是一项专门为诱导多能干细胞设计的、高灵敏度的细胞遗传学质控服务，是评估iPSC遗传稳定性与临床安全性的核心技术。

- 采用适用于iPSC的温和制片技术，并结合高分辨率G显带，获得清晰可辨的染色体带型
- 高内涵图像采集与分析：在全自动显微成像系统下分析至少20个分裂相，确保统计严谨性
- 专业报告与解读：出具详尽的核型分析报告，并对任何异常提供专业解读和建议

下单方式

关注道鹏生物微信公众号（xpbiomed）→进入右下角“检测服务”→立即下单

2025 Dr.Cell 技术服务列表

细胞鉴定检测项目	货号
Mycoplasma 支原体污染检测	T80100001
TR - 人源细胞STR鉴定	T80200001
细胞物种鉴定	T80300001
Sterility 无菌检测	T80400001
Cell Senescence 细胞衰老检测	T80500001
微生物测序（黑胶虫）	T80600001
核型分析 - 人源细胞	T80700001

外泌体检测项目		货号
NTA-光散射模式	粒径和浓度	T81010001
	电位和浓度	T81010010
	粒径、电位和浓度	T81010003
	粒径和浓度（NS300）	T81010004
纳米流式-粒径	粒径和浓度	T81010011
纳米流式-蛋白	纳米流式蛋白指标	T81010014
负染透射电镜（TEM）		T81010002
外泌体WB-1个蛋白指标		T81010007
外泌体WB-3个蛋白指标		T81010008
外泌体提取服务-超离		T81010006
外泌体提取服务-SEC		T81010013
BCA蛋白检测		T81010015
qPCR		T81010016
miRNA测序		T81010021
LncRNA测序		T81010022
蛋白组学（TMT）		T81010032
蛋白组学（4D label-free）		T81010034

企业概述

上海道鹏生物科技有限公司成立于2005年,由留美归国学者与来自斯坦福大学、罗氏等机构的顶尖专家共同创立,在细胞培养领域拥有超过20年的专业经验。公司秉承“始于科学,忠于科学”的理念,专注于科研试剂、耗材及相关服务的开发与推广。

2012年,创立自主品牌VivaCell,聚焦常规细胞培养试剂的生产与研发,产品广泛应用于细胞治疗、3D细胞培养(类器官)、外泌体、细胞遗传学等领域。同时,旗下Dr.Cell(希尔博士)还提供细胞检测服务-NTA, STR, Species ID, Senescence, Sterility等。

公司早在2016年就通过ISO13485医疗器械质量管理体系认证,凭借稳定可靠的产品质量,在羊水培养基、去外泌体胎牛血清、外泌体检测等细分领域取得全国领先的市场地位。未来,道鹏生物将持续致力于打造国际化本土标杆企业,助力中国生物医药产业实现进口替代与技术自强。

20 年+
专注细胞培养领域

18 国+
畅销全球18个国家

8200 平+
高标准生产基地

10000 个+
服务客户数量



ISO13485医疗器械质量体系认证



生产质量体系符合cGMP



欧盟CE认证 遵循国际标准

始于科学 忠于科学

始于科学 忠于科学

上海道鹏生物科技有限公司
Shanghai XP Biomed Ltd.

- 🏠 公司地址：上海市奉贤区汇丰西路2082号9号楼
- ☎ 咨询热线：400-820-3979
- ✉ 公司邮箱：info@xpbiomed.com
- 🌐 官方网址：www.xpbiomed.com



道鹏生物 官方微信

※ 请联系当地授权经销商购买正规产品